PolyBranch HCR™ miRNA、circRNA FISH Kit(Plant)

Catalog Number:HKR13Z


【Product Information】

Product Name	Catalog No	20T	Storage	Shelf Life
Wash Buffer (powder, dissolve in 10L ddH ₂ O)	HKR13-1	10L	Short-term:4℃ Long-term:-20℃	12 months
Proteinase K (100x)	HKR13-2	20uL		
Pectinase/Cellulase Mixture	HKR13-3	2mL		
AP-Mouse Anti-Digoxin (100x)	HKR13-4	20uL		
NBT Chromogen Solution A (20x)	HKR13-5	100uL		
NBT Chromogen Solution B (20x)	HKR13-6	100uL		
Blocking Buffer	HKR13-7	2mL		
Probe:	HKR13-8	2mL		
Pre-Amplification Probe	HKR13-9	2mL		
Amplification Probe	HKR13-10	2mL		

【Product Description】

The miRNA hybridization kit includes a combination probe, a pre amplification structure, and an amplification probe. When two probes bind to adjacent positions of the target sequence, a small segment of the top sequence can bind to the pre amplification structure. The repeated

sequence on the pre amplification structure can trigger a branched HCR to form a large nucleic acid aggregate. Each nucleic acid strand is labeled with digoxin, and the antibody recognizes digoxin for color development.

【Probe Information】

【Tissue Fixation】

Tissue Type	Animal	Plant	Frozen Samples	Cell Climb Slides	Cells
Treatment	Fix at RT for 12h, paraffin embed	Vacuum fix for 1h, RT fix for 12h, paraffin embed.	Dehydrate in 15% sucrose at 4°C for 8h, then in 30% sucrose at 4°C for 8h, OCT embed.	Fix at 4°C for 2h.	scrape off cells, fix in 4% PFA at 4°C for 2h, wash with PBS, agarose embed.
Type	mRNA	lncRNA	circRNA	miRNA	rRNA
Treatment	Fix at RT for 12h	Fix at RT for 12h (<300bp: 24h)	Fix at RT for 12h	Fix at RT for 12h	Fix at RT for 12h (<300bp: 24h)

【Storage and Shipping】

Ship on wet ice; store at -20 °C for long-term or at 4 °C for short-term use. Shelf life: 6 months.

【Protocol】

1. Deparaffinization and Rehydration

Immerse slides sequentially in: Xylene I (15 min) → Xylene II (15 min) → Xylene III (15 min) → 100% ethanol (10 min) → 90% ethanol (10 min) → 80% ethanol (10 min) → 70% ethanol (10 min) → Rinse with distilled water.

2. Enzyme Repair

After slides are completely dry, draw a hydrophobic circle around the tissue using a histology pen (recommended: HKR14P In Situ Hybridization Pen). Place slides horizontally in a hybridization oven or humidified chamber. Add 100 µL of Pectinase/Cellulase Mixture onto the tissue and incubate at 37 °C for 30 min. Rinse with distilled water to stop the reaction. Remove excess liquid, then add 100 µL of Proteinase K repair solution (1X) and incubate at 37 °C for 30 min. Rinse with distilled water to stop the reaction.

3. Blocking

Remove excess liquid from slides. Add 100 µL of Blocking Buffer per slide and incubate at 37 °C for 30 min. Wash once for 5 min. **Wash steps: Place slide racks in a wash tank, add wash buffer (ensure samples are submerged), and shake at 60 rpm for 5 min.**

4. Probe Hybridization

Remove excess liquid from slides. Add 100 µL of probe per slide and incubate at 37 °C for 3 h or overnight (maintain humidity to prevent drying). Wash 5 times, 5 min each. **Wash steps: As described above.**

5. Pre-Amplification Hybridization

Remove excess liquid from slides. Add 100 µL of Pre-Amplification Probe to each slide

and incubate at 37 °C for 3 h or overnight (ensure humidity). Wash 5 times, 5 min each. **Washing steps: As described above.**

6. Probe Amplification

Remove excess liquid from slides. Add 100 µL of Amplification Probe per slide and incubate at 37 °C for 1.5 h (maintain humidity). Wash 5 times, 5 min each. **Wash steps: As described above.**

7. AP-Mouse Anti-Digoxin

Remove excess liquid from slides. Add 100 µL of AP-Mouse Anti-Digoxin (1X) per slide and incubate at 37 °C in a humidified chamber for 40 min. Wash 5 times, 5 min each. **Wash steps: As described above.**


8. Chromogenic Reaction


Prepare the chromogenic solution by adding 50 µL of NBT Chromogen Solution A and 50 µL of NBT Chromogen Solution B to 1 mL of deionized water. Add 100 µL of the mixture per slide and incubate at RT for 30 min. Rinse with distilled water to stop the reaction. Air-dry slides and mount with mounting medium.

【Precautions】

1. For research use only.
2. Wear lab coats and disposable gloves for safety.

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