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Paraffin section immunohistochemical experiment report

1. Experimental principle

Immunohistochemistry is the application of the basic principle of immunology - antigen-antibody reaction, that is, the principle of specific combination of antigen and antibody, through chemical reactions to label antibody chromogenic agents (fluorescein, enzymes, metal ions, isotopes) to determine the antigen (polypeptide and protein) in tissue cells, and to locate, qualitative and quantitative research on it. These are called immunohistochemistry (IHC) or immunocytochemistry (ICC).

2. Laboratory equipment and reagents

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3. Experimental procedure

3.1 Dewaxing and hydration: Put the sections into xylene I 12 min - xylene II 12 min - anhydrous ethanol I 6 min - 95% Ethyl alcohol for 6 min - 85% Ethyl alcohol for 6 min, rinsing with distilled water.

3.2 Repair: The tissue sections were placed in a repair box filled with EDTA repair solution (PH9.0) and the antigen repair was carried out in a microwave oven. The repair procedure was: medium fire for 8min, ceasefire for 8min, and medium-low fire for 7min. During this process, excessive evaporation of buffer solution should be prevented and the pieces should not be dried. After natural cooling, the slide was placed in PBS (PH7.4) and washed by shaking on the decolorizing shaker for 3 times, 5min each time.

3.3 Block endogenous peroxidase : The sections were incubated with 3% hydrogen peroxide at room temperature for 25min. The slides were placed in PBS (PH7.4) and washed three times on a decolorizing shaker for 5min each time.

3.4 Draw a circle around the tissue: Use a histochemical pen to draw a small circle 3-4 mm apart from the tissue along the outer outline of the tissue, and then add enough PBS to ensure that the subsequent addition of blocking serum, primary antibody, secondary antibody, and color development agent can completely cover the tissue without flowing along the sections.

3.5 Serum blocking: The tissue was uniformly covered with 3%BSA in the tissue chemical circle and closed at room temperature for 30min.

3.6 Add primary antibody: Gently shake off the sealing solution, add PBS to the section in a certain proportion of primary antibody, and the section is placed flat in a wet box at 4°C for overnight incubation. (Add a small amount of water to the wet box to prevent the antibody from evaporating)

3.7 Adding secondary antibody: The slide was placed in PBS (PH7.4) and washed by shaking on the decolorizing shaker for 3 times, 5min each time. After the slices were slightly dried, the tissue was covered with the secondary antibody (HRP label) of the corresponding species of the primary antibody, and incubated at room temperature for 50min.

3.8 DAB color development: The slides were placed in PBS (PH7.4) and washed by shaking on the decolorizing shaker for 3 times, 5min each time. After the sections were slightly dried, the freshly prepared DAB color developing solution (prepared by diluent and concentrated solution 1000:50) was added into the circle. The color developing time was controlled under the microscope. The positive color was brown and yellow, and the section was washed with pure water to terminate the color development.

3.9 Restaining of nucleus: Hematoxylin restained the nuclei for 2-3min. Wash with tap water, hematoxylin differentiate solution differentiation for 2s, rinse with tap water, hematoxylin bluing solution back blue 15-30s, water rinse.

3.10 Dehydration and sealing: Put the sections into 75% Ethyl alcohol for 4 min - 85% Ethyl alcohol for 4 min - 95% Ethyl alcohol for 4 min - Anhydrous ethanol for 4 min - xylene I for 4 min - xylene II for 4 min, after slightly drying, , take the slices out of xylene for a little dry, and seal them with neutral gum.

3.11 Microscopy: The results are interpreted under a white light microscope.

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4. Interpretation of results

The nucleus of hematoxylin stain is blue, and the positive signal of DAB is brown-yellow.

5. Precautions

5.1 Pay attention to whether the slice dewaxing is thorough;

5.2 Don't let the slices dry out during the experiment;

5.3 Be careful of the injured tissue on the pipette suction head during experimental operation;

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Biotechnology Co., Ltd.